



Steel Lab Cabinets

Genie Scientific provides a complete line of powder coated steel lab cabinets coming in standard sizes and configuration all in a color of your choice. Like all Genie casework they are constructed from 18-gauge steel and welded construction requiring little maintenance.

Base cabinets form the foundation of most workspaces and workbenches in laboratories. This is why they need to withstand harsh chemicals and support the loads of heavy equipment. Powder coating not only gives you a smooth finish that is aesthetically pleasing but also a needed resistant layer against chemicals.

Genie Steel Casework Line Includes:

- Base cabinets
- Wall Cabinets
- Tall Cabinets
- Solvent Storage
- Corrosive Storage

Cabinet Features:

- Removable back panels
- Five-knuckle stainless steel door hinges
- Heavy Duty Full-Extension ball bearing slides
- Adjustable Shelves

- Powder coat finish in a color of your choice

View our casework catalog ([hyperlink to catalog](#))