

fume hood

walk-in
slimline model






The Genie Slimline floor mounted walk-in style fumehood comes in three standard sizes, and is perfect for rolling in large experiments.

- *Deal direct with the manufacturer
- *Chain driven - trouble free safety glass sash
- *Phenolic upgraded liner - standard
- *One on one technical support
- *Visit our facility
- *Ship crated nationwide
- *Knock down for shipping

fume hood

walk-in
slimline model

8' hood specifications shown
see below for additional hood details



MODEL #	EXTERIOR DIMENSIONS			INTERIOR DIMENSIONS			TWO TUBE LIGHT FIXTURE(S)	Shipping Weight
	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	SS COLLAR SIZE		
GE6WIFH	72"	40.5"	95"	63"	24"	12" dia.	one 4'	895
GE8WIFH	96"	40.5"	95"	87"	24"	two 10" dia.	two 3'	1185
GE10WIFH	120"	40.5"	95"	111"	24"	two 12" dia.	two 4'	1315

Custom Hood Sizes Available