

# GENIE

S C I E N T I F I C



CLASSIC LAB LINE  
STEEL CASEWORK & TABLES  
G2 SERIES

## where you'll find it

---

### **Casework**

Standing	2
Sitting	3
Upper	4
Tall	5

### **Welded Tables**

Frame Sizes & Part Numbering	6
------------------------------	---

### **Welded Workstations**

Frame Sizes & Part Numbering	7
------------------------------	---

### **Options and Accessories**

Counter Top Materials	8
Hanging Cabinets	8
Power Strips	8
Shelf Materials	9
Adjustable Height Legs	10
Casters	10

### **G2 Welded Table Quote Sheet**

11

### **Colors**



12

### **G2 Tables & Workstations**


- Rugged fully welded construction to support heavy loads and ultimate rigidity
- Made from heavy gauge steel tube with a powder coat finish
- Weight capacity of 7,000LBS (fixed height units)
- Legs have leveling feet for easy micro adjustments
- Choose from several counter top option and accessories
- All exposed welded are ground smooth or superior finish & aesthetics

### STANDING HEIGHT FEATURES


- 35" h X 22" d
- Powder coated steel
- 18 gauge construction
- 5 knuckle stainless steel hinges
- Adjustable leveling feet
- Standard aluminum pull




B3518-000  
 B3524-000  
 B3530-000  
 B3536-000  
 B3542-000  
 B3548-000




B3530-200  
 B3536-200  
 B3542-200  
 B3548-200




B3530-20S  
 B3536-20S  
 B3542-20S  
 B3548-20S



B3558-20S



B3530-220  
 B3536-220  
 B3542-220  
 B3548-220



C3514-100



B3530-140B  
 B3536-140B  
 B3542-140B  
 B3548-140B



B3530-140A  
 B3536-140A  
 B3542-140A  
 B3548-140A



B3530-060  
 B3536-060  
 B3542-060  
 B3548-060




B3530-080  
 B3536-080  
 B3542-080  
 B3548-080



B3530-040  
 B3536-040  
 B3542-040  
 B3548-040

## SITTING HEIGHT FEATURES

- 29" h X 22" d
- Powder coated steel
- 18 gauge construction
- 5 knuckle stainless steel hinges
- Adjustable leveling feet
- Standard aluminum pull




B2918-000  
B2924-000  
B2930-000  
B2936-000  
B2942-000  
B2948-000



B2918-B11  
B2924-B11




B2918-100  
B2924-100



B2930-200  
B2936-200  
B2942-200  
B2948-200



B2918-110  
B2924-110



B2930-220  
B2936-220  
B2942-220  
B2948-220



B2918-030  
B2924-030



B2930-060  
B2936-060  
B2942-060  
B2948-060




B2930-030  
B2936-030  
B2942-030  
B2948-030



B2914-100

### UPPER CABINET FEATURES



- 30”h X 13”d
- Powder coated steel
- 18 gauge construction
- 5 knuckle stainless steel hinges
- 1/4” clear glass
- Standard aluminum pull



	W3018-000 W3024-000 W3030-000 W3036-000 W3042-000 W3048-000		W3030-2GSF W3036-2GSF W3042-2GSF W3048-2GSF
	W3018-1SH W3024-1SH		W3030-2GS W3036-2GS W3042-2GS W3048-2GS
	W3018-1GH W3024-1GH		W3030-2GH W3036-2GH W3042-2GH W3048-2GH
	W3030-2SH W3036-2SH W3042-2SH W3048-2SH		W3030-2SS W3036-2SS W3042-2SS W3048-2SS

## TALL CABINET FEATURES

- 84”h X 21”d
- Powder coated steel
- 18 gauge construction
- 5 knuckle stainless steel hinges
- Adjustable leveling feet
- 1/4” clear glass
- Standard aluminum pull



### WELDED TABLE


- Fully welded table frame made from steel tube
- Comes standard with cross bar at rear of table
- All exposed welds are ground smooth for superior finish & aesthetics

### TABLE SIZES

Dimensions refer to steel frame Height, Width and Depth.

Counter top dimensions are not included, but typical to have 1" overhang on all sides.

<b>FRAME HEIGHT = A DIM</b>
29", 32", 35"
<b>FRAME WIDTH = B DIM</b>
34", 46", 58", 70", 94", 118"
<b>FRAME DEPTH = C DIM</b>
22", 28", 34", 46"



### PART NUMBERING SYSTEM

The example below shows a G2 Welded Table Frame that is 35"h x 70"w x 28"d.

**WT357028**

Welded Table

Frame Width  
34", 46", 58",  
70", 94", 118"

Frame Height  
29", 32", 35"

Frame Depth  
22", 28", 34", 46"

## WELDED WORKSTATIONS


- Fully welded workstation made from steel tube
- Rear integral upright is 84”H, slotted for adjustable height shelving
- Each workstation comes standard with (2) tiers of shelving
- 94” workstations and larger require additional center upright
- Comes standard with welded upright at rear of table
- All exposed welds are ground smooth for superior finish & aesthetics

## TABLE SIZES

Dimensions refer to steel frame Height, Width and Depth.

Counter top dimensions are not included, but typical to have 1” overhang on front and no overhang on back.


UPRIGHT HEIGHT
84”
FRAME HEIGHT = A DIM
29”, 32”, 35”
FRAME WIDTH = B DIM
34”, 46”, 58”, 70”, 94”, 118”
FRAME DEPTH = C DIM
23”, 29”, 35”, 47”



## PART NUMBERING SYSTEM

The example below shows a G2 Welded Workstation Frame with 84” uprights, 35”h x 70”w x 29”d.

WW84357029





## COUNTER TOP MATERIALS

- With G2 tables it is typical for tables and workstations to have the same counter top as the rest of the lab.
- Choose from Phenolic, Epoxy Resin, Laminate, Stainless Steel, Butcher Block

## HANGING CABINETS

- Frames can support G2 hanging cabinets

## POWER STRIPS

- Extruded aluminum channels factory precut and wired to your specifications
- Select channel size, type and number of outlets
- Contact Genie for additional information



## SHELF MATERIAL OPTIONS

- Workstations come standard with two-tier shelving
- Standard shelf bracket is bookend style, either 12", 15", 18" deep
- Option for lower shelf either half depth or full depth

PART# & MATERIAL	THICKNESS	EDGE	NOTES
PRS = Phenolic Resin Shelf	3/4"	Square Edge	Chemical Resistant
BBS = Butcher Block Shelf	1 3/4"	Square Edge	Oil Rub Finish
ESDS = ESD Laminate Shelf	1"	Square	Particle Board Core
LAMS = Standard Laminate Shelf	1"	PVC	Particle Board Core
PCS = Powder Coated Steel Shelf	1" 18ga	Square Edge	Core Matches Frame
SSS = Stainless Steel Shelf	1" 18ga	Square Edge	T304 #4 Finish



## ADJUSTABLE HEIGHT LEGS

- Allows for 6" of manual height adjustment (requires socket wrench)
- Available with casters

**Part #**  
**ADJ30-36**



## CASTERS

- Plate caster with 1200lbs capacity (per set)
- All four casters are swivel and locking

**Part #**  
**CASTERS**






**G2 Welded Table Frame:**

Frame Material: 2"x2" 14ga Steel Tube

Counter Top:

Frame color:


QTY	"A" TOP WIDTH (INCHES)	"B" TOP DEPTH (INCHES)	"C" TOP THICKNESS (INCHES)	"D" FRAME WIDTH (INCHES)	"E" FRAME DEPTH (INCHES)	"F" FRAME HEIGHT (INCHES)	"G" OVERALL HEIGHT (INCHES)





## MADE IN THE USA

WHITE WHITE  
wh121



ANTIQUÉ WHITE  
wh16



BEIGE  
bg02



SILK GRAY  
ral7044



LIGHT GRAY  
gr230




LIGHT GUNMETAL  
gr10



GRAY BLUE  
ral5008




PASTEL BLUE  
5024



BRILLIANT BLUE  
5007




PURPLE RED  
3004



TEXTURE SILVER  
gr673



SEMI GLOSS BLACK  
bk05



\*The standard colors above are a representation only; color may vary slightly from finished product.  
Powder coated steel samples are available upon request.

Genie Scientific, Inc.

☎ **800.545.8816** p 714.545.1838 f 714.641.0496

17430 Mt. Cliffwood Circle, Fountain Valley, CA 92708

www.geniescientific.com

MADE IN THE USA

